

xeo+

Designed for more

MULTI-APPLICATION AESTHETICS
IN A SINGLE POWERFUL PLATFORM

AMI technologies

Model, not actual patient.

Reach for more in your practice

WHERE PRECISION MEETS PASSION

Born from a 20+ year legacy, xeo+ combines versatility, ease-of-use, and cutting-edge features to meet the continuously evolving demands of a thriving practice.

With an Nd:YAG laser handpiece, multiple light-based handpieces, and numerous modalities to choose from, you can customize your options based on your clients' needs and your evolving requirements. This flexibility allows you to expand your service offerings now and in the future.

XEO+ FEATURES AND BENEFITS

MULTI-APPLICATION TECHNOLOGY

Treat more than 25 skin concerns with a single device. Streamline practice operations without the need for multiple device.[1]

OPTIMIZED ERGONOMICS

Redesigned handpieces for precision, speed, accuracy and ease-of-use during treatments.

COOLVIEW® HANDPIECE INNOVATION

The sapphire window provides continuous contact cooling for fast and safe treatments to enhance patient experience.

INTUITIVE GUIDED NAVIGATION

Modern interface with plug-and-play navigation and treatment settings to minimize learning curves and maximize productivity.

VERSATILE HANDPIECE OPTIONS

Attract a wide range of skin types and concerns by offering tailored treatments.



Expand your aesthetic horizons

xeo+ offers a wide range of handpieces to address multiple aesthetic and dermatological concerns. Choose to customize your device or build up on it as your practice grows and evolves.

TREATMENT OPTIONS

Options	Skin Type	Hair	Pigmentation	Vascular	Revitalization	Wrinkles	Scar Revision	Skin Tone	Skin Texture	Redness
CoolView®	I-VI	✓		✓*	✓	✓	✓	✓	✓	✓
Limelight®	I-IV		✓	✓	✓			✓		✓
Prowave	I-IV	✓								
Acutip™	I-IV		✓	✓						

* CoolView Vascular settings are only indicated for skin types I-IV.



Nd:YAG

COOLVIEW (1064 nm)


Treat a broad range of conditions with the wavelength pioneered by Cutera. Our legacy first-to-market Nd:YAG laser is optimized for hair removal, vascular lesions and Cutera's signature Laser Genesis™ skin renewal treatment.

- Large spot sizes up to 16 mm
- Fast treatment times
- Signature Laser Genesis Treatment

Signature Treatment

LASER GENESIS

Cutera's signature Laser Genesis treatment is safe for all skin types. This no downtime, quick laser facial helps to diffuse redness and improve the look of texture and scars while advancing overall skin quality.



Model, not actual patient.

INTENSE PULSED LIGHT (IPL) (500 -1100 nm)


xeo+'s comprehensive options of effective IPL handpieces provide a complete range of treatments addressing vascular and pigmented lesions, photoaging, and hair removal.



LIMELIGHT (520 - 1100 nm)

This signature IPL treats a variety of skin conditions including benign pigmented lesions, vascular lesions, and photo damage.


- Safe for skin types I-IV
- Age and sun spots
- Rosacea
- Broken capillaries
- Facial spider veins



PROWAVE (680 - 1100 nm)

Prowave is a cutting-edge, non-invasive light therapy that provides long-term or permanent hair reduction for all hair types, including finer hair on the face and body.

- Safe for skin types I-IV
- Stable long term or permanent hair reduction
- Large spot size for quick treatments
- Cooled tip for increased comfort



ACUTIP (500 - 635 nm)

Our exclusive Acutip IPL tool is designed for precise targeting of pigment, age and sun spots, discrete pigmented and vascular lesions, as well as small facial veins with a focus on precision and efficacy.

- Safe for skin types I-IV
- Small tip allows for ease of use in treatment area
- Refined approach to smaller treatment areas for safety
- Spot treatment capability

Designed to deliver results

Before | After 2 Treatments



COOLVIEW®

Photos courtesy of J. Lezaic, DO

Before | After 2 Treatments



LIMELIGHT®

Photos courtesy of T. Raymond, PA-C

Before | After 6 Treatments



LASER GENESIS™

Photos courtesy of N. Desai, DO


Before | After 6 Treatments



LASER GENESIS™

Photos courtesy of Vein Institute and Medical Spa, Inc.

Before | After 8 Treatments



PROWAVE

Photos courtesy of Vein Institute and Medical Spa, Inc.

Before | After 3 Treatments



LIMELIGHT® + LASER GENESIS™

Photos courtesy of N. Desai, DO

Before | After 3 Treatments



LIMELIGHT® + ACUTIP™

Photos courtesy of M. Kaga, MD

Before | After 3 Treatments



LIMELIGHT® + ACUTIP™

Photos courtesy of M. Kaga, MD

Individual Results May Vary.

TECHNICAL SPECIFICATIONS

Specs	CoolView	Limelight	Prowave	Acutip
Wavelength	1064 nm	520 - 1100 nm*	680 - 1100 nm*	500 - 635 nm*
Fluence Range	3 - 300 J/cm^2	<30 J/cm^2	5 - 35 J/cm^2	≤ 24 J/cm^2
Pulse Duration	0.1 - 300 ms	15 – 61 ms*	3.5 - 90 ms*	16 ms
Repetition Rate	≤ 10 Hz	≤ 2 Hz	≤ 2 Hz	≤ 1 Hz
Spot-size	3 -16 mm CIRCULAR SPOT	10 x 30 mm	10 x 30 mm	6.35 mm CIRCULAR SPOT
Cooling	Sapphire Contact Cooling			
Console Dimensions	14" W x 19" D x 34" H			
Console Weight	150 lbs.			
Console Electrical	100 -120 VAC/20A or 200 - 240 VAC/16A			

*Automatically set by Program Mode

References: 1 xeo+ Series Operator Manual. Cutera, Inc. D002931, REV. A, 05/24


Indications for Use *1064 nm Nd:YAG Indications:* The 1064 nm wavelength is indicated for the removal of unwanted hair; stable long-term, or permanent hair reduction through selective targeting of melanin in hair follicles; treatment of pseudofolliculitis barbae (PFB); coagulation and hemostasis of benign vascular lesions such as, but not limited to, hemangiomas, port wine stains, telangiectasias, rosacea, venous lakes, leg veins, spider veins and poikiloderma of Civatte; treatment of wrinkles such as, but not limited to, periocular and perioral wrinkles; reduction of red pigmentation in hypertrophic scars where vascularity is an integral part of the scar; treatment of warts, treatment of mild to moderate inflammatory acne vulgaris. The 1064 nm wavelength indicated for use on all skin types (Fitzpatrick I-VI), including tanned skin. *Optional Pulsed Light (IPL) Handpiece Indications:* ACUTIP and LIMELIGHT: The Acutip and Limelight handpieces are indicated for the treatment of vascular lesions, and treatment of benign pigmented lesions. PROWAVE: The Prowave handpiece is indicated for the removal of unwanted hair from skin types I-IV and stable long-term or permanent hair reduction. **Important Safety Information:** The following treatment-related expected transient side effects and possible adverse events may occur during or following treatment with xeo+: Erythema; Edema; Purpura; Hyperpigmentation; Hypopigmentation; Burns; Erosion; Epidermal crusting; Blistering; Scarring; Deep tissue injury; Prolonged wound healing; Temporary or permanent gray hair; Stimulation of terminal hair growth; Undesired hair loss in hair-bearing areas; Red rash/bumps; Hemosiderin staining; textural changes; Cutaneous indentations; Sun exposure, tanning beds, and artificial tanning may increase the risk of side effects and adverse events. Do not use xeo+ if you are pregnant or are undergoing treatment for skin cancer. Additional warnings and precautions can be discussed during a treatment consultation. **Contraindications:** The xeo+ series of systems is contraindicated for: pregnant patients; patients undergoing treatment for skin cancer. **Warnings:** Do not treat over dysplastic nevi or questionable pigmented lesions. Do not treat over open wounds. Do not treat over or close to tattoos or permanent make-up. Do not treat over or close to metal or electronic implants. Hair removal by lasers or intense pulse light sources can cause increased hair growth in some individuals. Based upon currently available data, the highest risk groups for this response are females of Mediterranean, Middle Eastern, and South Asian heritage treated on the face and neck. **Precautions:** Use caution when treating patients with any of the following: History of coagulopathies; History of keloids or hypertrophic scarring; History of herpes - pre-treatment with an antiviral may be indicated History of vitiligo - may cause de-pigmentation; Diabetes - may impede wound healing Metal or electronic implants; Current or recent use of the following medications: Accutane – do not treat if taken in the past 6 months; Gold Therapy – may cause blue-gray discoloration; Anticoagulants – may increase risk of purpura or bruising; Photosensitizing medications/agents. In addition, observe the following precautions: Use caution when treating over hair bearing areas. Laser energy may affect hair growth. Place moist gauze between the lips and teeth if treating near the mouth. Laser energy may affect teeth. Treat at least at least 15 cm (6 in) from pacemakers/ defibrillators. Reaction to metal or electronic implants is unknown. Refer to the Operator Manual to view the laser safety labels of the product.

Intended for North American audiences.

©2025 Cutera, Inc. All rights reserved. Model; Not actual patient. AP004117 rB (05/25)

About Cutera
Cutera was founded in 1998 by veteran laser and optical engineers who envisioned that technology should and could do more. We have since grown into a global leader of dermatology and aesthetics devices that appeal to forward-thinking clinicians who seek the next generation of performance, safety, and efficacy.

Cutera, Inc. World Headquarters | 3240 Bayshore Boulevard, Brisbane, California 94005 USA
888-4-CUTERA or 415-657-5500 | info@cutera.com | cutera.com



AMI technologies
ami.co.il | amiweb@ami.co.il | 09-7760300 משווק בלעדית ע"י אמי טכנולוגיות

מק"ט: CAT_XEO_P_0625_DR_EN | מהדורה: 0625